40 mm

FI

EOLAN

## ICHTHOPASTE™

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### Kuvaus

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Zinc Paste and Ichthammol Bandage

### Description

ICHTHOPASTE is a bandage made of open weave cloth impregnated with a paste containing 6.32% Zinc Oxide BP and 2% Ichthammol BP rolled onto a plastic core. It also contains Purified water, Glycerine, Gelatin (porcine), Emulsifying wax, Phenoxethanol. Each bandage is presented in a wrapper and sealed in a polythene bag contained in a cardboard carton.

### Mode of action

ICHTHOPASTE provides topical treatment. It provides a moist wound healing environment over a venous leg ulcer, helps to reduce skin irritation and soothes and protects surrounding skin. It will not dry out in use.

### Indications

ICHTHOPASTE is indicated to assist the management of venous leg ulcers. Where venous insufficiency exists, the paste bandage should be used under graduated compression bandaging. ICHTHOPASTE is also suitable for use in the management of chronic eczema/dermatitis, where occlusion is indicated.

### Contra-indications

Do not use in known cases of sensitivity or allergy to any of the ingredients.

### Precautions

Where venous insufficiency exists, ICHTHOPASTE should be used under graduated compression bandaging, after first assessing the patient to exclude arterial disease. The use of Doppler ultrasound is recommended for this purpose. Compression bandaging should not be used on patients who have arterial disease (ABI <0.8) and is not recommended for use on diabetic patients with advanced small vessel disease. Failure to detect reduced arterial flow can result in pressure necrosis, amputation, or even death.

### The risk of a patient having arterial as well as venous disease rises with age.

The skin of leg ulcer patients is easily sensitised to topical medicaments, including preservatives. Sensitisation should be suspected in patients particularly where there is deterioration of the surrounding skin. Such patients should be referred for specialist diagnosis, including patch testing. One of the functions of occlusive bandages is to increase absorption. Care should be taken, therefore, if it is decided to apply topical steroid, anaesthetic or retinoid preparations under the bandage, as their

### Recommended duration of use adults and children

Change bandage as required, can be worn up to 7 days.

### Possible side effects

Ichthopaste can cause side effects, although not everybody gets them. Has been reported:

- local skin reactions such as rash, reddening, itching and burning

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©Evolan Pharma AB
Trade Mark of Evolan Pharma AB
Evolan Pharma AB, Svärdvägen 19, Box 120, 182 12 Danderyd, Sverige
™Varemærke tilhørende Evolan Pharma AB
© Evolan Pharma AB